


Anti NF-L mAb 47:3 (UD1)

| | |
|--|--|
| Category | Mouse monoclonal |
| Purification | Protein G purified cell supernatant (serum free media). Purity ≥ 90 %. |
| Labelling | No label |
| Ig class | IgG1 _κ |
| Description/Specificity | UD1 recognizes the conserved rod domain of the neurofilament light protein (NF-L). |
| Application | ELISA SIMOA Immunohistochemistry Immunoblotting |
| Antigen Recognized in Species (tested so far) | Human, bovine, rat, mouse, sheep, macaque |
| Working Dilution | Application dependent |
| Storage Shelf-life | Long term at - 20°C. Shelf-life 1 year from shipping when stored at - 20°C. |
| Amount | 100 µg in PBS pH 7.4 |
| Article No. | 27016-100 |
| Origin | The product has been manufactured in Sweden. |



Address
UmanDiagnostics AB
Tvistevägen 48C
90736 Umeå
Sweden

Contact
Phone: +46 90 777 880
e-mail: info@umandiagnostics.se

Website
www.umandiagnostics.se

Issued: 2018-05