

Anti NF-L mAb 2:1 (UD2)

Category	Mouse monoclonal
Purification antibody	Protein G purified cell supernatant (serum free media). Purity ≥ 90 %.
Labelling	No label
Ig class	IgG1 _k
Description/Specificity	UD2 recognizes the conserved rod domain of the neurofilament light protein (NF-L).
Application	ELISA SIMOA Immunohistochemistry Immunoblotting
Antigen Recognized in Species (tested so far)	Human, bovine, rat, mouse, sheep, macaque
Working Dilution	Application dependent
Storage Shelf-life	Long term at - 20°C. Shelf-life 1 year from shipping when stored at - 20°C.
Amount	100 µg in PBS pH 7.4
Article No.	27017-100
Origin	The product has been manufactured in Sweden.



Address
UmanDiagnostics AB
Tvistevägen 48C
90736 Umeå
Sweden

Contact
Phone: +46 90 777 880
e-mail: info@umandiagnostics.se

Website
www.umandiagnostics.se

Issued: 2018-05