


## Biotin Labelled Anti NF-L mAb 2:1 (UD3)

Category	Mouse monoclonal
Purification antibody	Protein G purified cell supernatant (serum free media). Purity ≥90 %.
Labelling	Biotin
Ig class	IgG1 <sub>κ</sub>
Description/Specificity	UD3 recognizes the conserved rod domain of the neurofilament light protein (NF-L).
Application	ELISA SIMOA Immunohistochemistry Immunoblotting
Antigen Recognized in Species (tested so far)	Human, bovine, rat, mouse, sheep, macaque
Working Dilution	Application dependent
Storage Shelf-life	Long term at - 20°C. Shelf-life 1 year from shipping when stored at - 20°C.
Amount	100 µg in PBS pH 7.2 (10 mg/ml BSA & 0.01 % Sodium Azide)
Article No.	27018-100
Origin	The product has been manufactured in Sweden.



*Address*  
UmanDiagnostics AB  
Tvistevägen 48C  
90736 Umeå  
Sweden

*Contact*  
Phone: +46 90 777 880  
e-mail: [info@umandiagnostics.se](mailto:info@umandiagnostics.se)

*Website*  
[www.umandiagnostics.se](http://www.umandiagnostics.se)

Issued: 2018-09